

Galvani Systems MRI Eligibility Form

At the Time of the MRI Appointment:

- Review MRI scan conditions found at <https://galvani.bio/mri>
- Refer to the information below for MRI Scan Eligibility

Patient Name _____ Date of Birth _____

Physician name and office _____

Address, phone number _____

Galvani System IPG Model Number: _____ Serial Number: _____

Galvani System Lead Model Number: _____ Serial Number: _____

Date Eligibility was determined: _____

Patient and Galvani System Confirmation

- ☐ This patient does not have any additional Galvani components implanted other than the Galvani system being assessed for this MRI scan eligibility.




Has the Galvani System integrity been tested ?

- ☐ Galvani System Integrity confirms that the system is functioning nominally.
- ☐ Galvani System Integrity did not pass. MRI cannot not be performed.

Has Galvani System been programmed OFF ?

- ☐ The Galvani system stimulation is programmed OFF or set to Hibernation mode.

MRI Eligibility Mode

<input type="checkbox"/>	MR Conditional Head and Whole Body Scan Eligible (3 T)	
<input type="checkbox"/>	MR Conditional Head, Legs, and Feet Scan Eligible (1.5 T)	
<input type="checkbox"/>	The neurostimulation system MRI scan-type eligibility cannot be determined or is MR Unsafe	

Caution: Investigational device.
Limited by US Federal law to investigational use. Exclusively for Clinical Investigation (UK/EU).

Galvani Bioelectronics
Gunnels Wood Road
Stevenage, Hertfordshire, United Kingdom
Website: www.galvani.bio